



### Product description

|                     |  |
|---------------------|--|
| Order number        | 55.518   |
| Product description | Tube, working volume: 13 ml, (LxØ): 95 x 16.8 mm, material: PP, round base, transparent, push cap, without cap, 500 piece(s)/bag |

### Product characteristics

|                          |            |
|--------------------------|------------|
| Graduation               | no         |
| Closure type             | push cap   |
| Base shape               | round base |
| Volume of work           | 13 ml      |
| Centrifugation max (RCF) | 5000 x g   |

### Size

|                   |         |
|-------------------|---------|
| Diameter          | 16.8 mm |
| Length of product | 95 mm   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 13 ml, (LxØ): 95 x 16.8 mm, PP**

**Page 2**

**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

**Purity & certification**

|                           |                          |
|---------------------------|--------------------------|
| Satisfies the requirement | sealed at 0.6 bar/10 min |
| Product category          | no medical device   IVD  |
| Batched                   | no                       |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 36000         |
| Depth of case                 | 398 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 265 mm        |
| Case volume                   | 0.042 cbm     |
| Weight of product             | 0.0045 kg     |
| Weight of case                | 5.02 kg       |
| EAN of inner box              | 4038917082136 |
| EAN case                      | 4038917001229 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-28

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**