



### Product description

Order number	55.525.002
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, brown, suitable for Microvette® and micro tube with screw cap, 250 piece(s)/bag

### Product characteristics

Graduation	no
Closure type	push cap
Suitable for	Microvette® and micro tube with screw cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g

### Size

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product brown

### **Purity & certification**

Satisfies the requirement sealed at 0.6 bar/10 min  
Product category no medical device | IVD  
Batched no

### **Packaging**

Minimum order qty. 2000  
Type of inner packaging bag  
Piece(s) / inner box 250  
Piece(s) / case 2000  
Piece(s) / pallet 72000  
Depth of case 398 mm  
Width of case 398 mm  
Height of case 265 mm  
Case volume 0.042 cbm  
Weight of product 0.0025 kg  
Weight of case 5.57 kg  
EAN of inner box 4038917109147  
EAN case 4038917012355

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

