



Product description

Order number	55.525.061
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, transparent, push cap, with graduations, 100 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	13 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 5 ml, (LxØ): 75 x 13 mm, PP

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	4000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	4000
Piece(s) / pallet	48000
Case volume	0.1237 cbm
Weight of product	0.001 kg
Weight of case	5.85 kg
EAN of inner box	4038917276825
EAN case	4038917276818

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT