




Product Specification
Tube, 5 ml, (LxØ): 75 x 12 mm, PP

Page 1



Product description

Order number	55.526.305
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	12 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 5 ml, (LxØ): 75 x 12 mm, PP

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	4000
Piece(s) / pallet	64000
Case volume	0.0813 cbm
Weight of product	0.0024 kg
Weight of case	10.09 kg
EAN of inner box	4038917367318
EAN case	4038917367301

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

