



Product description

Order number	55.532.300
Product description	Tube, working volume: 4 ml, (LxØ): 75 x 11.5 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	4 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	11.5 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Tube, 4 ml, (LxØ): 75 x 11.5 mm, PP

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	80000
Case volume	0.0382 cbm
Weight of product	0.002 kg
Weight of case	4.77 kg
EAN of inner box	4038917367349
EAN case	4038917367332

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT