



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 55.535.300   |
| Product description | Tube, working volume: 3.5 ml, (LxØ): 55 x 12 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag |

### **Product characteristics**

|                          |            |
|--------------------------|------------|
| Graduation               | no         |
| Closure type             | push cap   |
| Base shape               | round base |
| Volume of work           | 3,5 ml     |
| Centrifugation max (RCF) | 5000 x g   |
| Autoclavable (max. °C)   | 121 °C     |

### **Size**

|                   |       |
|-------------------|-------|
| Diameter          | 12 mm |
| Height of product | 55 mm |
| Length of product | 55 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 3.5 ml, (LxØ): 55 x 12 mm, PP**

**Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | transparent        |

**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 1000          |
| Piece (s) / outer case        | 5000          |
| Piece(s) / pallet             | 140000        |
| Case volume                   | 0.0649 cbm    |
| Weight of product             | 0.0017 kg     |
| Weight of case                | 8.84 kg       |
| EAN of inner box              | 4038917367370 |
| EAN case                      | 4038917367363 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

