



Product description

Order number	55.539
Product description	Tube, working volume: 11.5 ml, (LxØ): 100 x 15.7 mm, material: PP, round base, transparent, push cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	11,5 ml
Centrifugation max (RCF)	4000 x g

Size

Diameter	15.7 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	42000
Depth of case	398 mm
Width of case	398 mm
Height of case	250 mm
Case volume	0.0396 cbm
Weight of product	0.0046 kg
Weight of case	5.2 kg
EAN of inner box	4038917082327
EAN case	4038917281997

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

