



### **Product description**

Order number	55.598.006
Product description	Tube, working volume: 8.5 ml, (LxØ): 75 x 15.7 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	8,5 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Diameter	15.7 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 8.5 ml, (LxØ): 75 x 15.7 mm, PP**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	56000
Case volume	0.0649 cbm
Weight of product	0.0024 kg
Weight of case	5.9 kg
EAN of inner box	4038917363136
EAN case	4038917363129

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**