



### **Product description**

Order number	57.462
Product description	Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PS, conical base, transparent, push cap, with graduations, 500 piece(s)/bag

### **Product characteristics**

Graduation	yes
Closure type	push cap
Base shape	conical base
Volume of work	12 ml
Centrifugation max (RCF)	2500 x g

### **Size**

Diameter	17 mm
Length of product	110 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0063 kg
Weight of case	6.98 kg
EAN of inner box	4038917081702
EAN case	4038917000819

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

