



Product description

Order number 57.477.500
Product description Tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PP, conical base, transparent, suitable for Pharmacia/LKB Gamma Counter and Gamma Counter Cobra from Packard, 1,000 piece(s)/bag

Product characteristics

Graduation no
Closure type push cap
Suitable for Pharmacia/LKB Gamma Counter and Gamma Counter Cobra from Packard
Base shape conical base
Volume of work 4,5 ml
Centrifugation max (RCF) 5000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Diameter	12 mm
Length of product	75 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	4000
Type of inner packaging	bag
Piece(s) / inner box	1000
Piece(s) / case	4000
Piece(s) / pallet	96000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0025 kg
Weight of case	10.59 kg
EAN of inner box	4038917107778
EAN case	4038917011945

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

