



Product description

Order number	57.512.001
Product description	Tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PP, conical base, transparent, push cap, without cap, 500 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Base shape	conical base
Volume of work	4,5 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	12 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 50000
Type of smallest subpackaging StackPack
Piece(s) / inner box 500
Piece (s) / outer case 2000
Piece(s) / pallet 120000
Depth of box 387 mm
Width of box 77 mm
Height of box 176 mm
Depth of case 398 mm
Width of case 358 mm
Height of case 175 mm
Case volume 0.0249 cbm
Weight of product 0.0021 kg
Weight of case 4.88 kg
EAN of inner box 4038917108799
EAN case 4038917072489

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

