



### **Product description**

Order number	57.512.001
Product description	Tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PP, conical base, transparent, push cap, without cap, 500 piece(s)/StackPack

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	conical base
Volume of work	4,5 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Diameter	12 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

### **Packaging**

Minimum order qty. 50000  
Type of smallest subpackaging StackPack  
Piece(s) / inner box 500  
Piece (s) / outer case 2000  
Piece(s) / pallet 120000  
Depth of box 387 mm  
Width of box 77 mm  
Height of box 176 mm  
Depth of case 398 mm  
Width of case 358 mm  
Height of case 175 mm  
Case volume 0.0249 cbm  
Weight of product 0.0021 kg  
Weight of case 4.88 kg  
EAN of inner box 4038917108799  
EAN case 4038917072489

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

