



### **Product description**

Order number	57.512
Product description	Tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PP, conical base, transparent, push cap, without cap, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	conical base
Volume of work	4,5 ml
Centrifugation max (RCF)	5000 x g

### **Size**

Diameter	12 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	4000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	4000
Piece(s) / pallet	96000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0021 kg
Weight of case	9.26 kg
EAN of inner box	4038917082075
EAN case	4038917001168

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

