



### **Product description**

Order number	58.485.001
Product description	Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PS, flat base, transparent, push cap, natural, sterile, 100 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	standard cap
Closure type	push cap
Base shape	flat base
Volume of work	7 ml
Centrifugation max (RCF)	2500 x g

### **Size**

Diameter	16 mm
Length of product	50 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 7 ml, (LxØ): 50 x 16 mm, PS**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	natural

**Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	32000
Depth of case	648 mm
Width of case	420 mm
Height of case	145 mm
Case volume	0.0395 cbm
Weight of product	0.004 kg
Weight of case	4.66 kg
EAN of inner box	4038917108027
EAN case	4038917012126

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

