



Product Specification
Tube, 21 ml, (LxØ): 65 x 23.5 mm, PS

Page 1



Product description

Order number	58.489
Product description	Tube, working volume: 21 ml, (LxØ): 65 x 23.5 mm, material: PS, flat base, transparent, push cap, without cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	flat base
Volume of work	21 ml
Centrifugation max (RCF)	1500 x g

Size

Diameter	23.5 mm
Length of product	65 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0057 kg
Weight of case	6.33 kg
EAN of inner box	4038917081962
EAN case	4038917001076

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

