



### **Product description**

Order number	58.618.008
Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PS, skirted conical base, brown, push cap, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	skirted conical base
Volume of work	12 ml

### **Size**

Diameter	16.5 mm
Height of product	95 mm
Length of product	95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tube, 12 ml, (LxØ): 95 x 16.5 mm, PS**

**Page 2**

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	brown

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Case volume	0.0492 cbm
Weight of product	0.0056 kg
Weight of case	6.1 kg
EAN of inner box	4038917115148
EAN case	4038917032414

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**