



Product description

Order number	60.506.024
Product description	Tube, working volume: 10 ml, (LxØ): 79 x 16 mm, material: PP, round base, transparent, screw cap, without cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	10 ml
Centrifugation max (RCF)	5000 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	16 mm
Length of product	79 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 10 ml, (LxØ): 79 x 16 mm, PP

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	32000
Case volume	0.0813 cbm
Weight of product	0.0032 kg
Weight of case	7.39 kg
EAN of inner box	4038917108713
EAN case	4038917035873

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT