



### **Product description**

Order number	60.540.316
Product description	Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16 mm, material: PP, round base, transparent, screw cap, natural, cap assembled, sterile, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	13 ml
Centrifugation max (RCF)	4000 x g

### **Size**

Diameter	16 mm
Length of product	101 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

**Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Purity standard	sterile
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	28000
Case volume	0.029 cbm
Weight of product	0.0058 kg
Weight of case	3.43 kg
EAN of inner box	4038917375986
EAN case	4038917375979

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

