



### Product description

Order number	60.542.530
Product description	Screw cap tube, working volume: 8 ml, (LxØ): 57 x 16.5 mm, material: PP, flat base, transparent, screw cap, without cap, 1,000 piece(s)/bag

### Product characteristics

Graduation	no
Closure type	screw cap
Base shape	flat base
Volume of work	8 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

### Size

Diameter	16.5 mm
Height of product	57 mm
Length of product	57 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	60000
Case volume	0.0282 cbm
Weight of product	0.0032 kg
Weight of case	3.67 kg
EAN of inner box	4038917111614
EAN case	4038917013260

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

