



Product description

Order number	60.548.012
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, transparent, HD-PE, without cap, with graduations, 250 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	50 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	28 mm
Length of product	114 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)

Purity & certification

Satisfies the requirement	ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	10000
Depth of case	598 mm
Width of case	398 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP

Page 2

Height of case	350 mm
Case volume	0.0833 cbm
Weight of product	0.0098 kg
Weight of case	5.72 kg
EAN of inner box	4038917111973
EAN case	4038917057868

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT