



Product description

| | |
|---------------------|--|
| Order number | 60.549 |
| Product description | Screw cap tube, working volume: 3.5 ml, (LxØ): 66 x 11.5 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap enclosed, 1,000 piece(s)/bag |

Product characteristics

| | |
|--------------------------|----------------------|
| Graduation | no |
| Closure type | screw cap |
| Base shape | skirted conical base |
| Volume of work | 3,5 ml |
| Centrifugation max (RCF) | 4000 x g |

Size

| | |
|-------------------|---------|
| Diameter | 11.5 mm |
| Length of product | 66 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | natural |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | IATA |
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 84000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 250 mm |
| Case volume | 0.0396 cbm |
| Weight of product | 0.0022 kg |
| Weight of case | 5.03 kg |
| EAN of inner box | 4038917082402 |
| EAN case | 4038917001502 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

