



Product description

Order number	60.550.009
Product description	Screw cap tube, working volume: 7 ml, (LxØ): 82 x 13 mm, material: PP, flat base, transparent, screw cap, natural, cap assembled, sterile, 100 piece(s)/bag

Product characteristics

Graduation	no
Closure type	screw cap
Base shape	flat base
Volume of work	7 ml
Centrifugation max (RCF)	4000 x g

Size

Diameter	13 mm
Length of product	82 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Purity standard	sterile
Batched	yes

Packaging

Minimum order qty.	50000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Case volume	0.0436 cbm
Weight of product	0.0044 kg
Weight of case	4.9 kg
EAN of inner box	4038917112246
EAN case	4038917032520

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

