



### Product description

Order number	60.551.001
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 79 x 16 mm, material: PP, round base with skirted base, transparent, screw cap, natural, cap assembled, sterile, 100 piece(s)/bag

### Product characteristics

Graduation	no
Cap	standard screw cap
Closure type	screw cap
Base shape	round base with skirted base
Volume of work	10 ml
Centrifugation max (RCF)	4000 x g

### Size

Diameter	16 mm
Length of product	79 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

**Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	413 mm
Width of case	257 mm
Height of case	464 mm
Case volume	0.0492 cbm
Weight of product	0.005 kg
Weight of case	5.6 kg
EAN of inner box	4038917112673
EAN case	4038917044103

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

