



### **Product description**

Order number	60.595.241
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, material: PP, skirted conical base, transparent, screw cap, white, cap assembled, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	standard screw cap
Closure type	screw cap
Base shape	skirted conical base
Volume of work	25 ml
Centrifugation max (RCF)	4500 x g

### **Size**

Diameter	25 mm
Length of product	90 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white

### **Purity & certification**

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	15000
Depth of case	458 mm
Width of case	408 mm
Height of case	260 mm
Case volume	0.0486 cbm
Weight of product	0.0086 kg
Weight of case	5.05 kg
EAN of inner box	4066916002473
EAN case	4066916002466

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

