



**Product Specification**  
**Screw cap tube, 8 ml, (LxØ): 94 x 14 mm, PC**

**Page 1**



**Product description**

Order number	60.9922.941
Product description	Screw cap tube, working volume: 8 ml, (LxØ): 94 x 14 mm, material: PC, round base, transparent, screw cap, yellow, cap assembled, sterile, 100 piece(s)/bag

**Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	8 ml
Centrifugation max (RCF)	8800 x g

**Size**

Diameter	14 mm
Length of product	94 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polycarbonate (PC)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	yellow

**Purity & certification**

Satisfies the requirement	ADR
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	608 mm
Width of case	324 mm
Height of case	264 mm
Case volume	0.052 cbm
Weight of product	0.0064 kg
Weight of case	7.2 kg
EAN of inner box	4038917153478
EAN case	4038917027847

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

