



Product description

Order number	62.470.300
Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PS, flat base, transparent, push cap, with print, label/print: white, with graduations, 500 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	flat base
Volume of work	12 ml

Size

Diameter	16.5 mm
Height of product	95 mm
Length of product	95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Case volume	0.0492 cbm
Weight of product	0.0056 kg
Weight of case	6.09 kg
EAN of inner box	4038917368001
EAN case	4038917367998

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

