



Product description

Order number	62.470
Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PS, flat base, transparent, push cap, without cap, with print, label/print: white, with graduations, 500 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	flat base
Volume of work	12 ml
Centrifugation max (RCF)	3500 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Diameter	16.5 mm
Length of product	95 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0055 kg
Weight of case	6.21 kg
EAN of inner box	4038917081788
EAN case	4038917000895

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

