



### **Product description**

Order number	62.476.020
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, with print, label/print: green, with graduations, 1,000 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Colour of print/label	green
Graduation	yes
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	3500 x g

### **Size**

Diameter	12 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	80000
Case volume	0.0382 cbm
Weight of product	0.0026 kg
Weight of case	6.2 kg
EAN of inner box	4038917107389
EAN case	4038917226073

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

