



Product description

Order number	62.486
Product description	Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PS, flat base, transparent, push cap, with print, label/print: white, with graduations, 1,000 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	flat base
Volume of work	7 ml
Centrifugation max (RCF)	2500 x g

Size

Diameter	16 mm
Length of product	50 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0028 kg
Weight of case	6.37 kg
EAN of inner box	4038917081931
EAN case	4038917001045

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

