



## Product description

Order number	62.503.210
Product description	Screw cap tube, working volume: 6.5 ml, Ø: 13 mm, round base, transparent, without cap, with print, label/print: white, with graduations, 500 piece(s)/bag

## Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Light protection	no
Base shape	round base
Volume of work	6,5 ml
Centrifugation max (RCF)	3000 x g

## Size

Diameter	13 mm
Height of product	90 mm

## Material & colours

Colour of product	transparent
Tube material	Polypropylene (PP)

## Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Safety data sheet	no
Batched	yes

## Packaging

Minimum order qty.	50000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	398 mm
Width of case	398 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 6.5 ml, Ø: 13 mm, with print**

Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0027 kg
Weight of case	6.12 kg
EAN of inner box	4038917429061
EAN case	4038917429054

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

