



### Product description

Order number	62.515.006
Product description	Tube, working volume: 13 ml, (LxØ): 100 x 16 mm, material: PP, round base, transparent, push cap, natural, cap assembled, with print, label/print: white, with graduations, sterile, 25 piece(s)/bag

### Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	Ventilation cap
Closure type	push cap
Base shape	round base
Volume of work	13 ml
Centrifugation max (RCF)	3000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

Diameter	16 mm
Length of product	100 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	natural

### Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	500
Piece(s) / pallet	24000
Depth of case	598 mm
Width of case	398 mm
Height of case	145 mm
Case volume	0.0345 cbm
Weight of product	0.0049 kg
Weight of case	3.05 kg
EAN of inner box	4038917108836
EAN case	4038917012263

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

