



Product description

Order number	62.515.028
Product description	Tube, working volume: 13 ml, (LxØ): 100 x 16 mm, material: PP, round base, transparent, push cap, natural, cap assembled, with print, label/print: white, with graduations, sterile, 1 piece(s)/blister

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Cap	Ventilation cap
Closure type	push cap
Base shape	round base
Volume of work	13 ml
Centrifugation max (RCF)	3000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Diameter	16 mm
Length of product	100 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	13500
Depth of box	351 mm
Width of box	151 mm
Height of box	171 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.0049 kg
Weight of case	4.06 kg
EAN of inner box	4038917233118
EAN case	4038917233088

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

