



### **Product description**

Order number	62.536.009
Product description	Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PP, flat base, transparent, push cap, with print, label/print: white, with graduations, 1,000 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	flat base
Volume of work	7 ml
Autoclavable (max. °C)	121 °C

### **Size**

Diameter	16 mm
Height of product	50 mm
Length of product	50 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Case volume	0.0492 cbm
Weight of product	0.002 kg
Weight of case	5.1 kg
EAN of inner box	4038917110037
EAN case	4038917033107

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

