



Product description

Order number	62.547.205
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, transparent, HD-PE, cap: red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 25 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	conical base
Volume of work	50 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 50 ml, (LxØ): 114 x 28 mm, PP,
with print**

Page 2

Size

Diameter 28 mm
Length of product 114 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap red

Purity & certification

Satisfies the requirement ADR
Product category not regulated
CE certified no CE mark
Purity standard free from DNA/DNase/RNase, free from
pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation Electron irradiation
Pyrogen-free yes
Non-cytotoxic yes
Batched yes

Packaging

Minimum order qty. 25000
Type of smallest bag
subpackaging
Piece(s) / inner box 25
Piece (s) / outer case 250
Piece(s) / pallet 9000
Depth of case 464 mm
Width of case 413 mm
Height of case 257 mm
Case volume 0.0492 cbm
Weight of product 0.0136 kg
Weight of case 4.15 kg
EAN of inner box 4038917111881
EAN case 4038917033275

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2026-03-10

This document was prepared by EDP support and is valid without signature.



SARSTEDT