



### **Product description**

Order number	62.547.254
Product description	Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, transparent, HD-PE, cap: red, cap assembled, with print, label/print: blue/white, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 25 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Colour of print/label	blue/white
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	conical base
Volume of work	50 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Diameter	28 mm
Length of product	114 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red

### **Purity & certification**

Satisfies the requirement	ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

### **Packaging**

Minimum order qty.	300
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	300
Piece(s) / pallet	8400
Depth of case	580 mm
Width of case	328 mm
Height of case	258 mm
Case volume	0.0491 cbm
Weight of product	0.0128 kg
Weight of case	4.61 kg
EAN of inner box	4038917111898
EAN case	4038917013420

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

