



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 62.547.354   |
| Product description | Screw cap tube, working volume: 50 ml, (LxØ): 114 x 28 mm, material: PP, conical base, brown, HD-PE, cap: red, cap assembled, with print, label/print: blue/white, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 25 piece(s)/bag |

### **Product characteristics**

|                          |              |
|--------------------------|--------------|
| Label/ Print             | with print   |
| Colour of print/label    | blue/white   |
| Graduation               | yes          |
| Closure type             | screw cap    |
| Base shape               | conical base |
| Volume of work           | 50 ml        |
| Centrifugation max (RCF) | 20000 x g    |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

|                   |        |
|-------------------|--------|
| Diameter          | 28 mm  |
| Length of product | 114 mm |

### Material & colours

|                   |                                   |
|-------------------|-----------------------------------|
| Product material  | Polypropylene (PP)                |
| Colour of product | brown                             |
| Closure material  | High Density Polyethylene (HD-PE) |
| Colour of cap     | red                               |

### Purity & certification

|                           |  |
|---------------------------|--|
| Satisfies the requirement | ADR  |
| Product category          | In vitro diagnostic device, CE   |
| CE certified              | CE   |
| Purity standard           | free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile |
| Sterilisation             | Electron irradiation   |
| Pyrogen-free              | yes  |
| Batched                   | yes  |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 300           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 25            |
| Piece (s) / outer case        | 300           |
| Piece(s) / pallet             | 8400          |
| Depth of case                 | 588 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 255 mm        |
| Case volume                   | 0.0597 cbm    |
| Weight of product             | 0.0128 kg     |
| Weight of case                | 4.54 kg       |
| EAN of inner box              | 4038917157759 |
| EAN case                      | 4038917157742 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

