



**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,**  
**with print**

**Page 1**



**Product description**

Order number	62.554.002
Product description	Screw cap tube, working volume: 15 ml, (LxØ): 120 x 17 mm, material: PP, conical base, transparent, screw cap, red, cap assembled, with print, label/print: white/blue, with graduations, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 50 piece(s)/polystyrene container

**Product characteristics**

Label/ Print	with print
Colour of print/label	white/blue
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	conical base
Volume of work	15 ml
Centrifugation max (RCF)	20000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 15 ml, (LxØ): 120 x 17 mm, PP,  
with print**

**Page 2**

**Size**

Diameter	17 mm
Length of product	120 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red

**Purity & certification**

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	polystyrene container
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	16000
Depth of case	435 mm
Width of case	143 mm
Height of case	702 mm
Case volume	0.0436 cbm
Weight of product	0.0057 kg
Weight of case	4.13 kg
EAN of inner box	4038917565103
EAN case	4038917565110

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

