



**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 57 x 15.3 mm, PP,  
with print**

**Page 1**



**Product description**

Order number	62.558.201
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 57 x 15.3 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap assembled, with print, label/print: white, Labelling area, with graduations, sterile, 100 piece(s)/bag

**Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	skirted conical base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Size**

Diameter	15.3 mm
Length of product	57 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

**Purity & certification**

Satisfies the requirement	IATA
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of case	598 mm
Width of case	398 mm
Height of case	145 mm
Case volume	0.0345 cbm
Weight of product	0.0037 kg
Weight of case	4.25 kg
EAN of inner box	4038917113472
EAN case	4038917280730

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

