



Product description

Order number	62.9924.272
Product description	Tube, working volume: 10 ml, (LxØ): 100 x 16 mm, material: PP, conical base, transparent, push cap, with print, label/print: white, with graduations, 100 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	conical base
Volume of work	10 ml
Centrifugation max (RCF)	4500 x g

Size

Diameter	16 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Depth of case	608 mm
Width of case	324 mm
Height of case	264 mm
Case volume	0.052 cbm
Weight of product	0.0042 kg
Weight of case	5 kg
EAN of inner box	4038917155120
EAN case	4038917037525

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

