



Product description

Order number	62.9924.282
Product description	Tube, working volume: 10 ml, (LxØ): 100 x 16 mm, material: PP, conical base, translucent, with printed graduations, label/print: white, with graduations, 1,000 piece(s)/bag

Product characteristics

Label/ Print	with printed graduations
Colour of print/label	white
Graduation	yes
Base shape	conical base
Volume of work	10 ml
Centrifugation max (RCF)	5000 x g
Centrifugation temperature	21
Centrifugation time	15 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 10 ml, (LxØ): 100 x 16 mm, PP, with printed graduations

Page 2

Size

Diameter	16 mm
Length of product	100 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	translucent
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	sealed at 18"hg/10min
Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	64000
Depth of case	548 mm
Width of case	373 mm
Height of case	408 mm
Case volume	0.0834 cbm
Weight of case	10 kg
EAN of inner box	4038917155175
EAN case	4038917398374

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-16

This document was prepared by EDP support and is valid without signature.



SARSTEDT