



## Product Specification

**Screw cap tube, 5 ml, 5 ml, (LxØ): 92 x 15.3 mm,  
conical false bottom, rounded tube bottom, PP,  
250 piece(s)/StackPack**

**Page 1**

### Product description

Order number	63.1592
Product description	Screw cap tube, working volume: 5 ml, nominal volume: 5 ml, (LxØ): 92 x 15.3 mm, conical false bottom, rounded tube bottom, transparent, material: PP, with paper thermo label, 250 piece(s)/StackPack, 1,000 piece(s)/carton

### Product characteristics

Label/ Print	with paper thermo label
Closure type	screw cap
Base shape	conical false bottom, rounded tube bottom
Volume of work	5 ml

### Size

Sample volume	5 ml
Diameter	15.3 mm
Length of product	92 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Tube material	Polypropylene (PP)

### Purity & certification

Product category	no medical device   IVD
Batched	yes

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	StackPack
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of box	336 mm
Width of box	95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 5 ml, 5 ml, (LxØ): 92 x 15.3 mm,**  
**conical false bottom, rounded tube bottom, PP,**  
**250 piece(s)/StackPack**

Height of box	177 mm
Depth of case	358 mm
Width of case	348 mm
Height of case	207 mm
Case volume	0.0258 cbm
Weight of product	0.0042 kg
Weight of case	4.93 kg
EAN of inner box	4038917269452
EAN case	4038917269445

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

