




Product Specification
Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP,
with paper label

Page 1



Product description

Order number	63.9922.248
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, material: PP, skirted conical base, transparent, screw cap, yellow, cap assembled, with paper label, label/print: white/yellow, Patient information, with graduations, sterile, 1 piece(s)/blister

Product characteristics

Label/ Print	with paper label
Colour of print/label	white/yellow
Graduation	yes
Closure type	screw cap
Base shape	skirted conical base
Volume of work	25 ml
Centrifugation max (RCF)	4500 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP,
with paper label

Page 2

Size

Diameter	25 mm
Length of product	90 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

Purity & certification

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	300
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	300
Piece(s) / pallet	4800
Depth of case	548 mm
Width of case	408 mm
Height of case	373 mm
Case volume	0.0834 cbm
Weight of product	0.0086 kg
Weight of case	4.3 kg
EAN of inner box	4038917238649
EAN case	4038917238632

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

