



Product description

| | |
|---------------------|--|
| Order number | 63.9922.254 |
| Product description | Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, material: PP, skirted conical base, transparent, screw cap, yellow, cap assembled, with paper label, label/print: white/yellow, Patient information, with graduations, sterile, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|----------------------|
| Label/ Print | with paper label |
| Colour of print/label | white/yellow |
| Graduation | yes |
| Closure type | screw cap |
| Base shape | skirted conical base |
| Volume of work | 25 ml |
| Centrifugation max (RCF) | 4500 x g |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP,
with paper label

Size

| | |
|-------------------|-------|
| Diameter | 25 mm |
| Length of product | 90 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | yellow |

Purity & certification

| | |
|---------------------------|-------------------------------|
| Satisfies the requirement | IATA |
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Gamma irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 15000 |
| Depth of case | 608 mm |
| Width of case | 324 mm |
| Height of case | 264 mm |
| Case volume | 0.052 cbm |
| Weight of product | 0.0084 kg |
| Weight of case | 4.98 kg |
| EAN of inner box | 4038917151573 |
| EAN case | 4038917027786 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

