




## Product Specification

**Multi-purpose container, max. volume: 500 ml,  
(LxØ): 150 x 70 mm, graduated, PP, transparent**

**Page 1**



### Product description

Order number	75.9922.812
Product description	Multi-purpose container, max. volume: 500 ml, (LxØ): 150 x 70 mm, transparent, cap: natural, graduated, material: PP, screw cap, cap assembled, 40 piece(s)/bag

### Product characteristics

Graduation	yes
Closure type	screw cap
Base shape	flat base
Volume of work	500 ml

### Size

Sample volume	500 ml
Max. volume	500 ml
Diameter	70 mm
Height of product	150 mm
Length of product	150 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Multi-purpose container, max. volume: 500 ml,  
(LxØ): 150 x 70 mm, graduated, PP, transparent**

**Page 2**

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural
Tube material	Polypropylene (PP)

### Purity & certification

Product category	no medical device   IVD
Latex-free	yes
Batched	yes

### Packaging

Minimum order qty.	80
Type of smallest subpackaging	bag
Piece(s) / inner box	40
Piece (s) / outer case	80
Piece(s) / pallet	1280
Depth of case	548 mm
Width of case	373 mm
Height of case	408 mm
Case volume	0.1445 cbm
Weight of product	0.0625 kg
Weight of case	5.87 kg
EAN of inner box	4038917153324
EAN case	4038917350921

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**