



**Product Specification**  
**Multi-purpose container, 420 ml, (LxØ): 150 x 70 mm, graduated, PP, transparent**

**Page 1**



**Product description**

Order number	75.9922.813
Product description	Multi-purpose container, max. working volume: 420 ml, (LxØ): 150 x 70 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, sterile, 40 piece(s)/bag

**Product characteristics**

Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Base shape	flat base
Volume of work	420 ml

**Size**

Sample volume	420 ml
Diameter	70 mm
Height of product	150 mm
Length of product	150 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	no medical device   IVD
Latex-free	yes
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

**Packaging**

Minimum order qty.	80
Type of smallest subpackaging	bag
Piece(s) / inner box	40
Piece (s) / outer case	80
Piece(s) / pallet	1280
Depth of case	548 mm
Width of case	373 mm
Height of case	408 mm
Case volume	0.0834 cbm
Weight of product	0.0655 kg
Weight of case	6.02 kg
EAN of inner box	4038917153331
EAN case	4038917350938

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

