



### **Product description**

Order number	93.894
Product description	Rack, material: PC, transparent, format: 6 x 3, (LxWxH): 137 x 70 x 40 mm, for 18 collection tubes, suitable for screw cap micro tubes, 1 piece(s)/inner box

### **Product characteristics**

Suitable for	screw cap micro tubes
Storage capacity	18
Format	6 x 3
Autoclavable (max. °C)	121 °C

### **Size**

Width of product	70 mm
Height of product	40 mm
Length of product	137 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polycarbonate (PC)
Colour of product	transparent

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1
Type of inner packaging	inner box
Piece(s) / inner box	1
Piece(s) / case	36
Piece(s) / pallet	2376
Depth of box	72 mm
Width of box	42 mm
Height of box	138 mm
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0,0182 cbm
Weight of product	0,1069 kg
Weight of case	4,91 kg
EAN of inner box	4038917084505
EAN case	4038917003315

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

