



**Product Specification**  
**Push cap, natural, suitable for micro tube Ø 10.8 mm**

**Page 1**



**Product description**

Order number	65.697
Product description	Push cap, natural, suitable for micro tube Ø 10.8 mm, flat, 1,000 piece(s)/bag

**Product characteristics**

Cap	standard cap
Closure type	push cap
Suitable for	micro tube Ø 10.8 mm

**Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Push cap, natural, suitable for micro tube Ø 10.8 mm**

**Page 2**

**Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	10000
Piece(s) / pallet	660000
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0.0182 cbm
Weight of product	0.0003 kg
Weight of case	3.31 kg
EAN of inner box	4038917083492
EAN case	4038917002417

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**