



Product description

Order number	65.710
Product description	Screw cap, natural, suitable for screw cap micro tube 72.733.201, for volume limitation to 100 µl, 500 piece(s)/bag

Product characteristics

Cap	flat cap
Closure type	screw cap
Suitable for	screw cap micro tube 72.733.201

Material & colours

Product material	Polypropylene (PP)
Colour of product	natural

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	220000
Depth of box	285 mm
Width of box	250 mm
Height of box	45 mm
Depth of case	568 mm
Width of case	348 mm
Height of case	155 mm
Case volume	0.0306 cbm
Weight of product	0.0006 kg
Weight of case	3.44 kg
EAN of inner box	4038917083607
EAN case	4038917080293

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

