



Product description

Order number	65.712
Product description	Screw cap, natural, suitable for screw cap micro tubes, 500 piece(s)/bag

Product characteristics

Cap	retaining loop cap
Closure type	screw cap
Suitable for	screw cap micro tubes

Material & colours

Product material	Polypropylene (PP)
Colour of product	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap, natural, suitable for screw cap micro tubes

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	330000
Depth of case	438 mm
Width of case	268 mm
Height of case	155 mm
Case volume	0.0182 cbm
Weight of product	0.0004 kg
Weight of case	2.51 kg
EAN of inner box	4038917083614
EAN case	4038917002516

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

