



**Product description**

|                     |  |
|---------------------|--|
| Order number        | 65.712   |
| Product description | Screw cap, natural, suitable for screw cap micro tubes, 500 piece(s)/bag |

**Product characteristics**

|              |                       |
|--------------|-----------------------|
| Cap          | retaining loop cap    |
| Closure type | screw cap             |
| Suitable for | screw cap micro tubes |

**Material & colours**

|                   |  |
|-------------------|--|
| Product material  | Polypropylene (PP)                     |
| Colour of product | natural                                |
| O-ring material   | Ethylene propylene diene rubber (EPDM) |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap, natural, suitable for screw cap micro tubes**

**Page 2**

**Purity & certification**

|                  |                                |
|------------------|--------------------------------|
| Product category | In vitro diagnostic device, CE |
| CE certified     | CE                             |
| Batched          | yes                            |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 5000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 5000          |
| Piece(s) / pallet             | 330000        |
| Depth of case                 | 438 mm        |
| Width of case                 | 268 mm        |
| Height of case                | 155 mm        |
| Case volume                   | 0.0182 cbm    |
| Weight of product             | 0.0004 kg     |
| Weight of case                | 2.51 kg       |
| EAN of inner box              | 4038917083614 |
| EAN case                      | 4038917002516 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-08

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**